



Latex Gloves

Supermax

General Consumables

Lightly powdered

Code	Size	Unit	Price/Unit	Price/Unit 10+	Price/Unit 30+
GLL1001	Small	Pk/100	\$7.20	\$6.85	\$6.45
GLL1002	Medium	Pk/100	\$7.20	\$6.85	\$6.45
GLL1003	Large	Pk/100	\$7.20	\$6.85	\$6.45
GLL1004	X-Large	Pk/100	\$7.20	\$6.85	\$6.45

Latex Gloves



Powder free - Microtouch

Code	Size	Unit	Price/Unit	Price/Unit 10+	Price/Unit 30+
GLLPF1001MT	Small	Pk/100	\$8.95	\$8.55	\$8.15
GLLPF1002MT	Medium	Pk/100	\$8.95	\$8.55	\$8.15
GLLPF1003MT	Large	Pk/100	\$8.95	\$8.55	\$8.15
GLLPF1004MT	X-Large	Pk/100	\$8.95	\$8.55	\$8.15

Nitrile Gloves

Supermax

General Consumables

Powder free - Microtouch

Code	Size	Unit	Price/Unit	Price/Unit 10+	Price/Unit 30+
GNBPF1001MT	Small	Pk/100	\$17.10	\$15.50	\$13.90
GNBPF1002MT	Medium	Pk/100	\$17.10	\$15.50	\$13.90
GNBPF1003MT	Large	Pk/100	\$17.10	\$15.50	\$13.90
GNBPF1004MT	X-Large	Pk/100	\$17.10	\$15.50	\$13.90

Nitrile Gloves



F10 Super Concentrate Disinfectant

F10

General Consumables

F10 Super Concentrate Disinfectant is a total spectrum disinfectant used for disinfection of hard surfaces, equipment and air spaces in pharmaceutical manufacturing, clean room and laboratory markets.

Code	Description	Unit	Price/unit
111	200ml F10SC Disinfectant	Each	\$23.50
113	1L F10SC Disinfectant	Each	\$93.00
103	5L F10SC Disinfectant	Each	\$398.00



F10 Disinfectant

pH Indicator Paper - Strips and Paper, Non-Bleeding

Macherey-Nagel

General Consumables

Code	pH	Description	Unit	Price/Unit	Price/Unit 10+
2401-100	pH 0-14	Indicator Strips	Pk/100	\$19.95	\$18.95
92115	pH 0-6	Indicator Strips	Pk/100	\$17.90	\$15.75
90203*	pH 1-11	Indicator Strips	Pk/100	\$17.90	\$15.75
90204*	pH 1-14	Indicator Paper	Pk/100	\$17.90	\$15.75
90205*	pH 0.5-5.5	Indicator Paper	Pk/100	\$17.90	\$15.75
90212*	pH 8-10	Indicator Strips	Pk/100	\$17.90	\$15.75
92118*	pH 2.0-9.0	Indicator Strips	Pk/100	\$17.90	\$15.75
92120*	pH 4.5-10	Indicator Strips	Pk/100	\$17.90	\$15.75
92122*	pH 6.0-10	Indicator Strips	Pk/100	\$17.90	\$15.75
92140*	pH 5.1-7.2	Indicator Strips	Pk/100	\$17.90	\$15.75



*Pricing strictly only applies to current stock. No backorders will be accepted or supplied at the price advertised.



Safety Eyewear - Eye Protection

Bolle

General Consumables

The Bolle safety eyewear glasses are designed in a sleek modern style and are ultra lightweight. All protection glasses are tested and approved to AS/NZS 1337 for medium impact protection.

Code	Model	Description	Unit	Price/Unit
1667201	Bandido	Clear frame/clear lens	Each	\$10.00
1668201	Blade	Clear frame/clear lens	Each	\$10.00
1671100	Mini B	Clear frame/clear lens	Each	\$13.00
1615501	Sidewinder	Gun metal frame/clear lens	Each	\$21.00
1661100	Extreme	Silver frame/clear lens	Each	\$26.50
1670102	Edge	Clear frame/contrast lens	Each	\$19.00
1614401	Prism	Clear frame/clear lens	Each	\$11.50
1650515	Override	Clear frame/clear lens	Each	\$10.00
1680303	Calibre	Clear frame/clear lens	Each	\$17.50




Bandido




Blade



Edge



Extreme



Sidewinder



Mini B

For the complete Bolle product range download the Bolle Flyer from the Lomb website.

Transfer Pipettes - Sterile & Non-Sterile

Samco

General Consumables

Samco transfer pipettes products are widely recognised for quality, safety and craftsmanship. Pipettes are available in bulk packs, individually sterile or sterile in packs of 20.

Sterile

Code	Description	Sterile	Unit	Pk/Ctn	Price/Unit	Price Ea/Ctn
222-1S	1ml	Individually wrapped	Pk/500	8	\$52.00	\$46.80
222-20S	1ml	Packs of 20	Pk/500	8	\$33.00	\$29.70
225-1S	7.5ml Graduated	Individually wrapped	Pk/500	8	\$51.00	\$45.90
232-20S	5.8ml Fine Tip	Packs of 20	Pk/500	8	\$39.15	\$37.20
251-1S	4.6ml SEDI-PET®	Individually wrapped	Pk/500	8	\$73.00	\$65.70
251-20S	4.6ml SEDI-PET®	Packs of 20	Pk/500	8	\$27.50	\$24.75
262-1S	6ml Extra Long 9"	Individually wrapped	Pk/400	8	\$106.50	\$95.85
262-20S	6ml Extra Long 9"	Packs of 20	Pk/400	8	\$45.00	\$40.50
282-20S	1.5ml Paediatric	Packs of 20	Pk/500	8	\$48.50	\$43.65



Samco Transfer Pipettes

Non sterile

222	1ml Graduated Large Bulb	Pk/500	10	\$17.00	\$15.30
225	3ml Graduated Large Bulb	Pk/500	10	\$17.00	\$15.30
232	5.8ml Fine Tip	Pk/500	10	\$33.00	\$29.70
234	8.7ml Fine Tip	Pk/400	10	\$33.00	\$29.70

Narrow Stem SEDI-PET®

241	1.2ml Short Tem	Pk/500	8	\$23.50	\$21.15
242	4ml Large Bulb	Pk/500	10	\$31.50	\$28.35



Foam Test Tube Racks

Bioplas

General Consumables

- Flame Retardant
- Charcoal colour
- Chemical and solvent resistant
- Float on H₂O bath with full tube and sample
- All tubes can be decanted simultaneously
- Low profile- will not tip over easily
- Lightweight- weighs less than two ounces

- Holds and locks 50 culture tubes securely
- Holds and locks 30 microcentrifuge and micro blood collecting tubes securely in racks
- More uniform incubation- tubes may be pushed down into H₂O bath while still being held securely in rack
- May be used in blood drawing trays-holds blood drawing tubes securely



Foam Test Tube Racks

Code	Tube Size Accommodated	Material	Capacity	Unit	Price/Unit
0010	10x75mm, 12x75mm & 13x100mm Tubes	Polyester	50 tubes/rack	Pk/6	\$50.40
0011	15x85mm, 16x100mm & 16x150mm Tubes	Polyester	50 tubes/rack	Pk/6	\$50.40
0020	250ul & 400ul Microcentrifuge Tubes	Polyester	30 tubes/rack	Pk/6	\$50.40
0021	500ul, 1.5ml & 1.8ml Microcentrifuge Tubes	Polyester	30 tubes/rack	Pk/6	\$50.40

Magnetic Stirring Bars - PTFE Coated

Cowie


General Consumables

Cylindrical

Code	Diam.	Length	Price/Ea
001.106	3mm	6mm	\$2.30
001.108	3mm	8mm	\$2.55
001.110.6	6mm	10mm	\$2.05
001.115.6	6mm	15mm	\$2.70
001.120	6mm	20mm	\$2.35
001.125.5	5mm	25mm	\$2.65
001.130	6mm	30mm	\$2.80
001.140	8mm	40mm	\$3.50
001.150	8mm	50mm	\$4.15
001.170	10mm	80mm	\$7.30

with pivot ring

Code	Diam.	Length	Price/Ea
001.313.8	8mm	13mm	\$2.55
001.315	4.5mm	15mm	\$2.30
001.315.8	8mm	15mm	\$2.65
001.320	6mm	20mm	\$2.60
001.325	6mm	25mm	\$2.75
001.330	6mm	30mm	\$3.15
001.340	8mm	40mm	\$4.15
001.350	8mm	50mm	\$4.85
001.370	10mm	70mm	\$ 6.45



Magnetic Stirring Bar

Lab Coats

Ivy Industries

General Consumables

Benefit and Feature Summary

- Stud front
- Three pockets
- Studded cuff
- Ply/cotton
- White.

Code	Size	Unit	Price/Unit
LABCOAT-WH (XS)	X-Small	Each	\$30.75
LABCOAT-WH(S)	Small	Each	\$30.75
LABCOAT-WH(M)	Medium	Each	\$30.75
LABCOAT-WH(L)	Large	Each	\$30.75
LABCOAT-WH(XL)	X-Large	Each	\$30.75



Lab Coats

Laboratory and Vacuum Tubing - Tygon® R-3603

Saint-Gobain

General Consumables

The most consistently reliable tubing for the transfer of liquids and gases, Tygon® Laboratory Tubing handles virtually all inorganic chemicals found in today's laboratories. It is non-oxidising and non-contaminating. Long-lasting and crack-resistant, Tygon® Laboratory Tubing is specially formulated to resist flex-fatigue and abrasion and is less permeable than rubber tubing. It features extra heavy walls that will withstand a full vacuum at room temperature (759mm of mercury at 23 °C and up to 686 mm of mercury at 60 °C). Coils are marked at 1-foot intervals for easy measuring.

Benefit and Feature Summary

- Outstanding chemical resistance
- Lot-to-lot consistency for reproducible results
- Smooth, polished inner wall
- Non-aging and non-oxidising



Discover more: View the product range and pricing by downloading this file at www.lomb.com.au/Tygon-R3603



Tygon Laboratory and Vacuum tubing is available in a wide range of sizes to suit virtually all possible applications.