



Cytofuge 12

The StatSpin CytoFuge 12 is a small, versatile and economically-priced cytocentrifuge designed to gently deposit cells in a monolayer on to microscope slides. The Cytofuge12 handles samples from routine body fluid preps to research methodologies involving application of multiple cell concentrations or controls on the same slide. A broad range of selectable speed and run times enables precise replication of most procedures developed on other cytocentrifuge systems.

Disposable Filter Concentrators

Save time and improve results by eliminating the need to clean and sterilise reusable sample chambers. **Sample volume range:** 50 to 400µL

Reusable/Disposable Cell Concentrators

Available in 1 and 3-well formats. These silicone gasket-sealed devices are made of transparent plastic and are designed for single use.

Sample volume range: 1-well cell concentrator 300 to 800µL
1-well cell concentrator 50 to 200µL

Benefit and Feature Summary

- Fast, Simple, Economical
- Prepares up to 12 Monolayer Cell Preps
- Suitable for use with a range of body fluids
- Excellent safety design
- Easy to program with up to 24 protocols
- Flexible speed, time and acceleration
- Choice of filter or gasket-seal concentrator formats

Technical data

Code	CF12*
Safety Features	- Cover interlock system - Leak-resistant rotor with transparent lid - Low voltage drive system (24 V) - Designed to meet requirements of 1EC-1010-2-020
Rotor Capacity	One to twelve slides
Speeds (Forces)	200-2,000 rpm 28-560 x g
Operating Times	1-99 minutes
Dimensions (H x W x D)	20.3 x 33.0 x 40.6cm
Weight	10kg

StatSpin

Cytology



Cytofuge 12

Code	CF12*
Description	StatSpin Cytofuge 12
Unit	Each
Price/Unit	\$8,553

Supplied complete with a 12 place rotor and start-up kit.

WEB

Discover more:

To download flyers, brochures and promotions simply visit: www.lomb.com.au/StatSpin

Fronine Cytology Stains

For over 25 years pathologists have relied on the consistent quality of Fronine brand stains and solvents in pathology laboratories throughout Australia. With a product range that extends to over 300 individual products, the Fronine brand includes stains and solvents for use in Histology, Microbiology, Haematology, Cytology and Parasitology.

Now you have the opportunity to purchase from a select range of routine Fronine stains and solvents at special promotional pricing and to also receive a FREE Cytology product with every purchase over \$250.

This is a strictly limited offer so you will need to act quickly to avoid missing out on the special promotional pricing and your FREE product.

Fronine

Cytology



Fronine Cytology Stains

Special pricing on Cytology stains

Code	Description	Size	Unit	Price/Unit
GG001	Fields Stain A	2.5L	Each	\$13.30
GG002		5L	Each	\$25.50
GG003	Field Stain B	2.5L	Each	\$13.30
GG004		5L	Each	\$25.50
II031	Papanicolaou EA36	2.5L	Each	\$61.20
II032		5L	Each	\$97.75
II037	Papanicolaou EA65	2.5L	Each	\$61.20
II038		5L	Each	\$97.75
II040	Papanicolaou OG6	2.5L	Each	\$61.20
II041		5L	Each	\$97.75

Bonus Offer

Purchase over \$250 worth of cytology stains from the list, at the advertised price, and receive* a pack of microtome blades for frozen tissue for FREE.

*Simply quote "FREE-CYTO" when placing your order.

PDF

Discover more: Download the complete Fronine Promotion at www.lomb.com.au/Fronine-Stains

*Non stock line. Please check delivery time when placing order.