




Test Tube Racks - Two/Three Tier

Kartell

Plasticware

Code	Description	Size	Unit	Price/Unit
KAR370	Two tier holds 10x18mm tubes	250x60x74mm	Each	\$ 10.45
KAR562	Three tier holds 24x20mm tubes	375x65x95mm	Each	\$ 18.95
KAR559	Three tier holds 12x12mm tubes	190x60x80mm	Each	\$ 11.50
KAR560	Three tier holds 12x20mm tubes	190x60x80mm	Each	\$ 11.50




Test Tube Racks - Universal, Stackable

Kartell

Plasticware

Code	Description	No. of holes	Colour	Unit	Price/Unit
KAR564-03	Test Tube Racks	90 x 13mm Ø	White	Each	\$ 13.50
KAR565-03	Test Tube Racks	60 x 16mm Ø	White	Each	\$ 13.50
KAR565-04	Test Tube Racks	60 x 16mm Ø	Blue	Each	\$ 13.50
KAR565-06	Test Tube Racks	60 x 16mm Ø	Yellow	Each	\$ 13.50
KAR566-03	Test Tube Racks	40 x 20mm Ø	White	Each	\$ 13.50
KAR566-10	Test Tube Racks	40 x 20mm Ø	Red	Each	\$ 13.50
KAR568-03	Test Tube Racks	24 x 30mm Ø	White	Each	\$ 14.90
KAR568-04	Test Tube Racks	24 x 30mm Ø	Blue	Each	\$ 14.90



Wide Neck Bottles - Graduated, Square

Kartell

Plasticware

Code	Volume	Grad.	Diameter (mm)	Mouth diam.	Unit	Price/Unit
KAR608	25ml	5	32 x 34 x 52mm	18mm	Each	\$1.35
KAR609	50ml	10	38 x 38 x 70mm	24mm	Each	\$1.70
KAR610	100ml	20	42 x 48 x 90mm	34.5mm	Each	\$2.05
KAR611	250ml	50	57 x 60 x 110mm	34.5mm	Each	\$2.80
KAR612	500ml	100	70 x 80 x 138mm	45mm	Each	\$3.40
KAR613	1ltr	100	80 x 103 x 176mm	58mm	Each	\$5.20



Eye Wash Station

Kartell

Plasticware

Eye wash bottle art. 383 to be ordered separately, is made of polyethylene provided with an eye cup with snap on cap, carefully designed to allow a quick and immediate flushing to remove foreign matters and prevent serious eye damage. The head has an air-vent valve to assure air refill of free space in the bottle thus preventing draw back of contaminated wash into the bottle. A perforated spray plate in the cup centre breaks up the stream to give the eye a gentle wash. A drain tube releases waste.

Code	Description	Unit	Price/Unit
KAR2386	Emergency Eye Wash Station	Each	\$41.85
KAR383	Integral Eye Wash Bottle M45	Each	\$18.90



Slide Staining System - Easy Dip™

Simport

Plasticware

The EasyDip™ Slide Staining Rack will hold up to 12 microslides, which fit into individual slots for free passage and rapid drainage of staining fluids. This system has two components: a square staining jar and a 12-position vertical slide rack. *EasyDip™ Slide Staining systems are available in **Green** (-12G), **Pink** (-12P), **White** (-12W) and **Yellow** (-12Y) and a Pk of assorted colours is also available.

Code	Description	Colour	Unit	Price/Unit
M900*	EasyDip™ Slide Staining Jars	Choose	Pk/6	\$99.50
M900-12AS#	EasyDip™ Slide Staining System	Assorted	Pk/6	\$129.50
M905-12DGY	EasyDip™ Slide Staining Rack	Dark Gray	Pk/6	\$103.95



* Six jars of the same colour are supplied without a staining rack. Staining racks must be purchased separately.
Each kit contains five jars (one of each colour) and one rack (M905-12DGY)



Weight Tray - High Impact, Graduated, Square, PS

Kartell

Plasticware

Code	Capacity	Dimensions	Unit	Price/Unit
KAR11103	7ml	43 x 43 x 8mm	Pk/1000	\$108.00
KAR11104	100ml	78.5 x 78.5 x 23mm	Pk/1000	\$207.00
KAR11105	250ml	134 x 134 x 23mm	Pk/1000	\$405.00



Weight Trays

Beakers - Low Form PP & PMP

Kartell

Plasticware

PP Moulded Graduations

Manufactured from polypropylene. Autoclavable to 121°C for 20min. Permanent graduations are moulded in for reproducibility. Excellent chemical resistance.

Material: Polypropylene



Beakers - PP Moulded

Code	Volume	Height	Unit	Price/Unit
KAR1801	25ml	49mm	Each	\$1.85
KAR1802	50ml	60mm	Each	\$2.05
KAR1803	100ml	72mm	Each	\$2.35
KAR1805	250ml	95mm	Each	\$2.90
KAR1806	500ml	119mm	Each	\$4.50
KAR1808	1ltr	147mm	Each	\$7.65
KAR1809	2ltr	183mm	Each	\$12.15
KAR811	3ltr	201mm	Each	\$23.40
KAR812	5ltr	229mm	Each	\$35.10

PP with Blue Graduations

Manufactured from polypropylene. Autoclavable to 121°C for 20min. Permanent graduations are moulded in for reproducibility. Excellent chemical resistance.

Material: Polypropylene



Beakers - with blue graduations

Code	Size	Height	Unit	Price/Unit
KAR1821	25ml	49mm	Each	\$2.15
KAR1822	50ml	60mm	Each	\$2.50
KAR1823	100ml	72mm	Each	\$2.70
KAR1824	250ml	95mm	Each	\$3.40
KAR1825	500ml	119mm	Each	\$5.40
KAR1826	1ltr	147mm	Each	\$9.90
KAR1827	2ltr	183mm	Each	\$13.95

Aspirator Bottles - Heavy Duty

Kartell

Plasticware

Designed for life time service. Can be used in sub zero temperatures. Very sturdy having a 2mm wall thickness whilst retaining translucency. The 3/4" BSP tap (supplied) mates with a 12 mm deep thread in the aspirator for total leakproof performance. Recommended that a few drops of vaseline be put on the thread before assembly.

Code	Volume	Mouth diam.	Diameter	Height	Handles	Unit	Price/Unit
KAR1660	5ltr	46mm	163mm	335mm	1	Each	\$36.00
KAR1662	10ltr	55mm	210mm	414mm	1	Each	\$45.00
KAR1664	25ltr	55mm	279mm	525mm	2	Each	\$82.80
KAR1666	50ltr	55mm	366mm	603mm	2	Each	\$127.80

Code	Description	Unit	Price
KAR375	Spigot 3/4" BSP thread	Each	\$17.10



Aspirator Bottles