



Compact Balances

***EK-i***  

---

***EW-i***  
Series



*EK-i*



*EW-i*



**AND** ...Clearly a Better Value  
A&D Company, Limited  
<http://www.aandd.jp>

# Choose the best for your weighing application

## High resolution up to 1/60,000! (EK-610i, EW-12Ki)

A wider selection of capacities & resolutions, and remote-zeroing has been included in a low-profile rugged design. For years the EK/EW series balances have been top selling products. Now the new EK-i Compact Balances continue their tradition with enhanced capacities and resolutions, and EW-i series comes with attractive triple range weighing. All feature RS-232C as standard and are designed to meet weights and measures regulations of various countries.



- Wide selection of 16 different models 30g ~ 12kg
- Large LCD display (16mm height) with back light – easy to read even at a wide viewing angle or dark environment
- Control zero, mode selection and data output with standard RS-232C
- Triple range weighing for EW-i series
- 9 built-in weighing units plus counting & percentage with user priority start-up
- Counting function with Automatic Counting Accuracy Improvement
- Standard comparator function – HI/OK/LO annunciator
- Full Digital Calibration with optional user definable mass setting
- Conformity to GLP with I.D. and serial number outputs
- Auto Power OFF
- Portable-compact size, light weight, rugged body
- Sealed key panel protection against dust and spills
- Low profile and small footprint for any space requirement
- Overload protection
- Stainless steel weighing pan
- Easy to install rechargeable Ni-MH battery (Optional)
- Optional underhook weighing for EK-4100i/6100i/6000i/12Ki, EW-12Ki
- Leveling feet and level bubble for accurate weighing
- Security ring to prevent theft

# EW-i



ions from a wider selection.

# EK-6100i)

and housing.  
 have been designed to continue  
 es of weighing.  
 untries (e.g. NTEP, EC).



# EK-i

## WEIGHING RANGE

(EW-i series only)  
 Triple range indication

## ACAI

Automatic Counting  
 Accuracy Improvement

## HI/OK/LO


Comparator Annunciator

## SAMPLE

Register a sample count or  
 100% in percentage

## PRINT

Transmit data to printer  
 or computer




### Multi weighing units

- |     |               |     |           |
|-----|---------------|-----|-----------|
| oz  | — Ounce (av)  | g   | — Gram    |
| ozt | — Troy ounce  | ct  | — Carat   |
| tl  | — Tael        | mom | — Momme   |
| lb  | — Pound (av)  | pcs | — Pieces  |
| GN  | — Grain unit  | %   | — Percent |
| dwt | — Pennyweight |     |           |

## GLP

Balance I.D. number allows  
 calibration output verification



## Options & Accessories

### Dot Matrix Compact Printer

\*AD-8121B



### Optional Ni-MH Pack

\*EKW-09i

For complete portability and use  
 in areas with power problems



### Carrying Case

\*EKW-12i



