


DURAN® Erlenmeyer Flasks with Baffles

Why use baffled flasks?

- Reduced cell growth due to limited oxygen transfer
- Baffled flasks
 - cause a turbulent flow
 - increase the gas exchange surface
 - produce a higher oxygen transfer

Benefits of DURAN® Erlenmeyer baffled flasks

- automated one-step production
 - = reproducible geometry
 - = comparable results
- complete with membrane screw cap



Order information:

Cat. No.	Description	Pack/Quantity
21 283 36 5	DURAN® baffled flask, GL 45, with 4 bottom baffles, complete with membrane screw cap, 250 ml	4
21 283 44 5	DURAN® baffled flask, GL 45, with 4 bottom baffles, complete with membrane screw cap, 500 ml	4

You can obtain these products from your specialist laboratory ware dealer.