



Life Science

Lab Water Systems

Cascada Water Purification Systems

Faster, consistent water purification with realtime quality monitoring for analytical applications

Designed for general laboratory and life sciences applications.
18.2 Megohm.cm (M½-cm) water for:

Applications

- ELISA tests
- Affinity chromatography
- Electrophoresis
- Media preparation
- Nucleic acid and protein isolation
- Monoclonal antibody production



- **Adjustable flow rate.**
- **Cartridge identity tags prevent errors**
electronically tagged cartridges prevent incorrect cartridge fitting and potential leaks. Tag data is displayed on the control monitor showing cartridge number, installation date, and recommended change out date. Information can be downloaded for validation.
- **No need for calibration modules**
unit calibrates electronic circuitry upon set up.
- **Comprehensive sanitisation**
unlike many competitive systems, Cascada's entire fluid pathway is included in the automated disinfection process ensuring complete sanitisation. Data tagged sanitisation blocks prevent other operations during sanitisation process.
- **Permanent record of data**
RS232 port allows for data collection and permanent record of water quality and system parameters—essential for compliance with Good Laboratory Practices (GLP) guidelines.
- **Fits anywhere**
compact design for benchtop, under bench or wall mounting.

Order Details

Code	Description	Unit
PAL-CAXXXROM2	Cascada RO system	Each
PAL-CAXXXBIOM2	Cascada Bio Polisher	Each
PAL-CAXXXANIM2	Cascada AN Polisher	Each
PAL-CAXXXLSM2	Cascada LS Polisher	Each
PAL-CAXXXLXM2	Cascada IX Polisher	Each
PAL-A653	35 Litre Docking Vessel	Each
PAL-A644	25 Litre Reservoir	Each
PAL-A645	40 Litre Reservoir	Each
PAL-A646	75 Litre Reservoir	Each

Contact your local Lomb Account Manager or Customer Service Representatives today to find out which Cascada System best suits your needs.

Call: 1300 725 973
Email: sales@lomb.com.au

• Chemicals • Consumables • Instrumentation



Lomb Scientific Pty Ltd

26 Atkinson Rd (PO Box 2223) Taren Point 2229 NSW

NSW 02-8536 6222 VIC 1300 725 973 QLD 07-3850 8466 WA: 08-9358 1844 All other states: 1300 725 973 Fax 1300 30 30 25