



EJ Series Compact Balance

A&D

Instrumentation

A&D's new compact balances, the EJ Series, offer an outstanding list of features, sure to provide you with excellent value for money.

Features include a unique and eye catching colour and stackable design facilitating easy storage. The units are also suitable for field work where battery operation and a high-visibility display are essential.

Benefit and Feature Summary

- Unique and eye-catching colour design
- 10 models available
- Your choice of multiple weighing units
- Counting and percentage functions included
- Anti-theft lock feature for added security
- Stackable body allows convenient storage for multiple units
- Calculates density automatically using the optional density kit
- Simple keyboard function
- Backlit display - ideal for poorly lit areas



EJ-1500 Square Pan Compact Balance



EJ-120 Round Pan Compact Balance

Code	EJ-120	EJ-200	EJ-300	EJ-410	EJ-610	EJ-1500	EJ-2000	EJ-3000	EJ-4100	EJ-61000
Weight capacity	120g	210g	310g	410g	610g	1500g	2100g	3100g	4100g	6100g
Resolution	0.01g					0.1g				
Weighting pan size	110mm Ø					127 x 140mm				
Price/Unit	\$495	\$540	\$590	\$640	\$690	\$495	\$540	\$590	\$640	\$690

WEB Discover more: Download flyers and promotions at www.lomb.com.au/AND

Omni TH_Q - Digital Tissue Homogeniser

Omni

Instrumentation

The new Omni TH_Q is a powerful tissue homogeniser designed to homogenise a wide range of samples from 0.03mL to 100mL. The TH_Q can usually process a tissue sample in less than 30 seconds and is compatible with a wide variety of generator probes for flexibility in application requirements.

The Omni TH_Q, Digital Tissue Homogeniser includes: TH_Q, motor, stand clamp, 6 soft tissue & 6 hard tissue Omni tip™ plastic homogeniser probes, tool kit, instruction manual and storage case.

Benefit Summary

- Integrated digital speed display & control for accurate, repeatable results (8,000 - 35,000 rpm)
- Whisper drive technology ensures ultra-quiet operation (<68db)
- Uses Omni Tip Plastic probes or Stainless Steel Probes
- Ergonomic handheld design or can be post mounted
- Energy efficient 20 watt motor
- 3 years motor warranty
- Processing range from 0.03mL-100mL
- Variable motor speeds from 5,000 to 35,000 rpm



PDF Discover more:
www.lomb.com.au/THQ



Omni TH_Q Homogeniser



Omni TH_Q Kit

Code	Description	Unit	Price/Unit
12-200-220*	Omni TH _Q Package	Each	\$2,942

Orbital Genie™ Benchtop Shaker

Scientific Industries

Instrumentation

Scientific Industries has introduced the new Orbital Genie™ Benchtop Orbital Shaker. The most versatile Orbital Genie™ is available from the makers of the Vortex-Genie®

Benefit and Feature Summary

- Choice of three different units depending on platform preference, ratcheting universal clamps or standard flask clamps or adhesive mat
- Variable shaking speed and option to reverse orbital action
- Programmable gentle ramp up to speed, so open containers will not spill
- Designed for use in cooled, standard and CO₂ incubators



Orbital Genie™

Code	Description	Unit	Price/Unit
SI-1506*	Orbital Genie™ Shaker, Ratcheting Clamps, 230V	Each	\$3,947
SI-1606*	Orbital Genie™ Shaker, Standard Flask Clamps, 230V	Each	\$3,023
SI-1706*	Orbital Genie™ Shaker, Adhering Mat, 230V	Each	\$3,287

PDF

Discover more:
www.lomb.com.au/OrbitalGenie

*Non stock line. Please check delivery time when placing order.



Low-Temperature Chillers - L-Series

PolyScience

Instrumentation

PolyScience L-series chillers are environment-friendly, economical alternatives to tap-water cooling. The unique design of PolyScience Compact Chillers maximises precious bench space without compromising cooling and pumping power.

Benefit and Feature Summary

- High performance at low temperatures
- Simple installation, operation, and maintenance
- Space-saving design
- User-settable high and low temperature and low flow rate alarms

An ideal and low cost companion for rotary evaporators



Chiller

Specifications

Model	LS5 Space Saving Chiller			LM6 Mini Chiller	
	M2 Centrifugal	MY Centrifugal	TY Turbine	GY Centrifugal	MY Centrifugal
Cooling capacity 20°C	1160	1190	830	520	390
	1030	1070	750	390	280
	680	700	500	250	170
	435	460	330	140	110
	230	240	140	-	-
Max. pump pressure	5.5psi 0.38bar	10.5psi 0.72bar	32psi 2.21bar	4.4psi 0.30bar	12.5psi 0.86bar
	47 x 24.4 x 59.4cm			48.9 x 24.4 x 41.9cm	
Overall dimensions					
Code	*LS52M21AA10E	*LS52MY1AA10E	*LS52TY1AA10E	*LM62GY1AA10E	*LM62MY1AA10E
Price	\$4,759	\$4,959	\$5,296	\$3,610	\$4,013

PDF

Discover more: Download the specifications at www.lomb.com.au/LChillers

Pipette Accuracy Testers

A&D

Instrumentation

The new range of Pipette Accuracy Testers from A&D provides everything you need to determine the inaccuracy and imprecision of any type or brand of pipette, routinely used in a laboratory.

Evaporation Trap

One of the challenges in weighing a small quantity of liquid (e.g. 50 µL or less) is controlling the environment to minimise errors due to loss of evaporation. The evaporation trap maintains high humidity inside and prevents evaporation of the test liquid. It is no longer necessary to set up and adjust the humidity of an entire room.

WinCT-Pipette

The special WinCT-Pipette Software provides complete data management and test results for all pipette calibrations. In accordance with quality assurance requirements, test results can be exported to an A4 or letter-size report format and printed.

Tester Kit Components

- Instruction manual
- Evaporation trap
- Liquid thermometer
- WinCT-Pipette (CD-ROM)
- Carrying case with a shoulder belt and a key
- Calibration weight with a pair of tweezers
- Sample cup with a holder (30 mL x 2 / 5 mL x 2)
- Balance incl. the weighing pan unit, breeze break, AC adaptor and AC adaptor ID label
- USB communications kit (USB converter, RS-232C cable, Instruction manual)

WEB

Discover more: Visit the A&D webpage and download the flyer at www.lomb.com.au/AND



Pipette not included!



Choose from three different testers

Specifications

Model	Weighting capacity	Repeatability	Unit	Price/Unit
*AD-4212B-PT	110g	0.1mg / 0.05mg	Each	\$8,995
*AD-4212A-PT	110g	0.15mg	Each	\$7,200
*FX-300i-PT	320g	1mg	Each	\$3,250



*Non stock line. Please check delivery time when placing order.